Hypertension

Human Fever

Asthma

Flexible Heating

**Europe / Middle-East / Africa**

Microlife AG  
 Espenstrasse 139  
 9443 Widnau / Switzerland  
 Tel. +41 / 71 727 70 30  
 Fax +41 / 71 727 70 39  
 Email admin@microlife.ch  
 www.microlife.com

**Asia**

Microlife Corporation.  
 9F, 431, RuiGang Road, NeiHu  
 Taipei, 114, Taiwan, R.O.C.  
 Tel. 886 2 8797-1288  
 Fax.886 2 8797-1283  
 Email service@microlife.com.tw  
 www.microlife.com

**North / Central / South America**

Microlife USA, Inc.  
 1617 Gulf to Bay Blvd., 2nd Floor Ste A  
 Clearwater, FL 33755 / USA  
 Tel. +1 727 442 5353  
 Fax +1 727 442 5377  
 Email msa@microlifeusa.com  
 www.microlife.com







IB WS 100 V16 3211

**Microlife WS 100**

- |                |                |
|----------------|----------------|
| <b>EN</b> → 2  | <b>TR</b> → 50 |
| <b>FR</b> → 6  | <b>GR</b> → 54 |
| <b>ES</b> → 10 | <b>AR</b> → 58 |
| <b>PT</b> → 14 | <b>FA</b> → 62 |
| <b>RU</b> → 18 |                |
| <b>UA</b> → 22 |                |
| <b>PL</b> → 26 |                |
| <b>HU</b> → 30 |                |
| <b>DE</b> → 34 |                |
| <b>SV</b> → 38 |                |
| <b>FI</b> → 42 |                |
| <b>NL</b> → 46 |                |

- ① SET Button
- ② User Button
- ③ Down Button
- ④ Up Button
- ⑤ Battery Compartment
- ⑥ KG / LB Switch

**Display**

- ⑦ Low Battery Indication
- ⑧ Zero Setting
- ⑨ Weight Reading
- ⑩ User Indication
- ⑪ Height
- ⑫ Age
- ⑬ Gender
- ⑭ Fitness Level
- ⑮ Body Fat
- ⑯ Body Fat Reading
- ⑰ Body Water
- ⑱ Body Water Reading
- ⑲ Error

Read the instructions carefully before using this device.

Dear Customer,

This new electronic personal scale allows you to evaluate your body composition and to automatically calculate the body fat and water percentage in function of your height, age and sex. The exact knowledge of your body composition (calculation of lean body mass, fat body mass and total body water contents), allows a correct estimate of the excess fat and at the same time determines the advisable weight as a function of your nutritional condition. Furthermore, the knowledge of your body composition allows to combine a correct nutrition and physical activity as essential conditions to improve the quality of life, promote health, and prevent pathologies. It has a memory function with a 10 person capacity. Capacity 150 kilos, divided into 100 grams.

Please read through these instructions carefully so that you understand all functions and safety information. We want you to be happy with your Microlife product. If you have any questions or problems, please contact Microlife-Customer Service. Your dealer or pharmacy will be able to give you the address of the Microlife dealer in your country. Alternatively, visit the Internet at [www.microlife.com](http://www.microlife.com) where you will find a wealth of invaluable information on our products.

Stay healthy – Microlife AG!

**Table of Contents**

- 1. Important Facts about Body Composition**
    - Body fat mass percentage
    - Water percentage
  - 2. Suggestions for Use**
    - Fitness
  - 3. Using the Instrument for the First Time**
    - How to insert the battery
    - KG / LB Switch
    - Use as simple scales
  - 4. Inserting Personal Data**
  - 5. Calculation of Weight, Body Fat and Water Percentage**
  - 6. Safety, Care, Accuracy Test, and Disposal**
    - Safety and protection
    - Instrument care
    - Disposal
  - 7. Guarantee**
  - 8. Technical Specifications**
- Guarantee Card (see Back Cover)**

Name of Purchaser / Nom de l'acheteur / Nombre del comprador / Nome do comprador / Ф.И.О. покупателя / П.І.П. покупця / Imię i nazwisko nabywcy / Vásárló neve / Name des Käufers / Inköparens namn / Ūstajan nimi / Alicinin Adı / Ονοματεπώνυμο αγοραστή / اسم المشتري / نام خریدار

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## 1. Important Facts about Body Composition

It is important to know the body composition of an individual in order to be able to estimate his/her health condition. It is not enough to know one's weight but it is necessary to define its quality, which means dividing the human body in its two major component masses:

- Lean body mass: bones, muscles, viscera, water, organs, blood
- Fat body mass: adipose cells

Under the same conditions of weight, these two masses have different volumes.

By means of electrodes placed on the platform, the scale sends a weak low intensity electric current that passes through the body mass. The electric current reacts in a different way according to the masses it passes through: the lean body mass opposes a weak resistance, while the fat body mass has a greater resistance. This body resistance, called bioelectrical impedance, varies also in function of sex, age and height.

For a correct control of one's ideal weight, it is important, besides the weight, to compare the body fat percentage calculated by the scale, with the parameters contained in the table and the body water percentage, with the information listed hereinafter.

**If the percentages are outside the normal values, we suggest you consult your physician.**

In order to better understand the value of your body fat percentage, we suggest you write down this value for a given period of time. In order to obtain a valid comparison, you should weigh yourself always at the same time of the day and under the same conditions.

### Body fat mass percentage

The percentage of fat listed in the following table represents the total fat, calculated by summing up the essential fat (or primary fat) and the spare fat (or storage fat). The essential fat is necessary since it is needed for the metabolism.

The percentage of this fat is different between the two sexes: about 4 % for men and about 12 % for women (percentage calculated on the total body weight).

The spare fat, which is theoretically superfluous, is necessary for those who carry out physical activities. Even in this case, the fat percentage is different between the two sexes: 12 % for men, 15 % for women (percentage calculated on the total body weight).

Therefore, the normal percentage of fat mass on the body weight is approximately 16 % for men and 27 % for women.

Such percentages differ with the person's age as shown in the table here below.

Women	Excellent	Good	Normal	Over-weight	Obese
Age	Data in %				
≤ 19	17.0	17.1-22.0	22.1-27.0	27.1-32.0	> 32.1
20-29	18.0	18.1-23.0	23.1-28.0	28.1-33.0	> 33.1
30-39	19.0	19.1-24.0	24.1-29.0	29.1-34.0	> 34.1
40-49	20.0	20.1-25.0	25.1-30.0	30.1-35.0	> 35.1
≥ 50	21.0	21.1-26.0	26.1-31.0	31.1-36.0	> 36.1

Men	Excellent	Good	Normal	Over-weight	Obese
Age	Data in %				
≤ 19	12.0	12.1-17.0	17.1-22.0	22.1-27.0	> 27.1
20-29	13.0	13.1-18.0	18.1-23.0	23.1-28.0	> 28.1
30-39	14.0	14.1-19.0	19.1-24.0	24.1-29.0	> 29.1
40-49	15.0	15.1-20.0	20.1-25.0	25.1-30.0	> 30.1
≥ 50	16.0	16.1-21.0	21.1-26.0	26.1-31.0	> 31.1

### Water percentage

Body water % is an important indication of the overall «well being» of a human being as water is the single most important component of body weight and represents more than half (about 60 %) of the total weight.

Research shows that 57 % body water is optimum, but this depends on the individual and normally there is no problem if the reading is 10 % higher or lower.

The body water % is greater in childhood (at birth it varies from 70-77 %) and it decreases with age (old people might have a body water of 45-55 %). Women (usually having more adipose tissue than men), have a water ratio of about 55-58 % while men have about 60-62 % of it.

With the increase of fat deposits, obese people also have a lower body water % than lean people. In general, body water % is an overall composite number. Below 40 % would normally be seen as too low, and above 70 % as too high. Human beings excrete a given quantity of water every day. This, however, must be regularly replaced, a water loss of about 10 % of the body weight, may put our health condition at risk.

So giving the body water % as a reading on the body fat scale gives an indication to the user if there are possible problems, in which case they would probably wish to consult their doctor.

Body weight variations that take place in a few hours are completely normal and are due to water retention, since the extra-cellular water (which includes interstitial fluid, plasma, lymph and transcellular fluid) is the only compartment subject to such rapid variations.

## 2. Suggestions for Use

- Place the scale on a flat, hard surface. Soft, inaccurate flooring (e.g. rugs, carpets, linoleum) can cause the scale to give uneven reading of your weight.
- Step onto the scale with your feet parallel and weight equally distributed.
- Stand still while the scale measures your weight.
- It is important to weigh yourself with dry, bare feet, correctly positioned on the electrodes. Moreover, for the scale to work properly, avoid your legs touching.
- The scale is safe to use as a simple scale by everyone being powered by a single 3 volt lithium battery. However, the following people should not perform a body fat or body water measurement:
  - pregnant women
  - people on cardiovascular medication
  - those fitted with a pace-maker or other forms of electronic, medical equipment
  - those with a temperature
  - children under the age of 10
  - those in dialysis, with oedemas on their legs, people with dimorphism or suffering from osteoporosis

- Always weigh yourself on the same scale each day at the same time, preferably undressed and before breakfast. To get the best results from your scale, weigh yourself twice, and if the two weights are different from each other, your weight is between the two readings. Wait for about 15 minutes after getting up, so that the water can distribute itself around your body.
- The surface of the scale can be slippery, if wet.
- If an error occurs during the analysis of your body fat and water percentage, the scale will turn off itself automatically.
- The precision of this fat and water percentage measurement can be altered in the following cases:
  - a high alcohol level
  - a high level of caffeine or drugs
  - after an intense sports activity
  - during an illness
  - during pregnancy
  - after heavy meals
  - in presence of dehydration problems
  - your personal data (height, age, sex) have not been input correctly
  - if you have wet or dirty feet

 **WARNING:** The body fat and water percentage is to be considered only as approximate information. For any further information, please consult your physician or dietician.

 **WARNING:** The measurement function of body fat and water percentage may be utilized only by people that belong to an age group that goes from 10 to 100 and that are between 100 and 200 cm tall.

### Fitness

In order to determine your body fat percentage in the most correct way, this scale will take into consideration not only your weight, height, age, and sex, but it will also take into account your lifestyle and the level of the physical activity that you practice.

For this reason, when you input your personal data you should choose among the 3 different FITNESS LEVELS.

#### • Fit 1:

This category includes all sedentary people who in the last months have not carried out any physical activity or have practised it at moderate level, less than twice a week for less than 15 minutes during each session.

This level may be selected by those people who want to start a diet together with physical activity aimed at reducing the fat body mass. The Fit 1 level is to be utilized for approximately 6/10 weeks from the start of your training program/diet.

#### • Fit 2:

This category includes all people who regularly practice physical activities (from 3 to 5 times per week, at least 20 minutes during each session). The Fit 2 level is to be utilized approximately for 8/12 weeks.

#### • Fit 3:

This category includes people (athletes, sportsmen) who carry out an intense daily physical activity (5 times or more per week) and that have developed flexibility, resistance, and strength in time.

## 3. Using the Instrument for the First Time

### How to insert the battery

This scale operates with one replaceable lithium battery. When the battery is dead the display shows «LO» ⑦ or the battery symbol. Insert a new battery in the battery compartment ⑤ with the positive pole upwards.

Do not discard the batteries together with the normal wastes, discard with special wastes.

Remove the battery if the scale is not used for long periods of time.

### KG / LB Switch

Found at the bottom of the scale, this switch ⑥ allows you to choose the unit of measure for weighing: kg, lb or st.

Select your desired unit of measure always when the scale is turned off (otherwise the unit of measure cannot be changed).

### Use as simple scale

1. Turn the scale on by pressing the centre of the platform.
2. Wait a few seconds and when «0.0» ⑧ appears on the display screen, the scale is ready for use.
3. Step onto the scale and wait without moving; after a few seconds, the display screen will show your weight ⑨.
 

Do not lean on anything while being weighed as this could alter the reading.
4. Step off the scale, the display screen will continue to show your weight for a few seconds then it will turn off automatically.

## 4. Entering Personal Data

This scale can memorise the data of 10 different persons.

1. Turn the scale on by pressing the centre of the platform.
2. Wait a few seconds and when «0.0» ⑧ appears on the display, press the SET button ①.
3. «P-01» ⑩ will flash on the display screen, indicating the «01» person: if you wish to change the memory location (memories from 01 to 10), press the «▼» ③ or «▲» ④ symbols. Once you have selected the memory location, press SET ① to confirm.
4. Insert your height (from 100 cm to 200 cm) ⑪ using the «▼» ③ or «▲» ④ symbols to increase or decrease the value. The display screen will give the value in cm or ft. Press SET ① to confirm.
5. Use the «▼» ③ or «▲» ④ symbols again to insert age (from 10 to 100) ⑫ and then press SET ① to confirm.
6. Continue by pressing the «▼» ③ or «▲» ④ symbols to define your sex ⑬. Press SET ① to confirm.
7. At this point, the display will show a fitness number «Fit1», «Fit2» or «Fit3» (depending on what fitness number was previously selected in the memory for the same user number) ⑭.
8. Now the scale has memorized all your data; the display will show them in sequence in order for you to verify that they have been input correctly. At this point the scale will turn off by itself.

To change the data entered, recall the relevant memory location and repeat the whole procedure.

To enter another person's data, select a free memory location and repeat the whole procedure described above.

 **WARNING:** When you are inside the function «personal data input» the scale will turn off automatically if no button is pressed within 10 seconds.

## 5. Calculation of Weight, Body Fat and Water Percentage

1. Turn the scale on by pressing the centre of the platform ②.
2. Wait a few seconds and when «**0.0**» ⑧ appears on the display, press the USER button ②. Press the USER button ② again until the memory location required appears. The data inserted will be viewed on the display screen and then «**0.0kg**» ⑧ will appear.
3. Step on the scale making sure you place your bare feet over the two electrodes. The display will show the following data in sequence:
  - weight ⑨
  - the «**bF %**» (bodyfat) sign ⑮
  - the percentage of fat in your body ⑯
  - the «**H2O %**» (water) sign ⑰
  - the percentage of water in your body ⑱Compare the percentage of bodyfat (bF) obtained with the parameters listed in the «Body fat mass percentage» table: if it falls outside the normal values, we suggest you consult your physician. Estimate the percentage of water (H<sub>2</sub>O) as displayed in the paragraph «Water percentage».
4. The scale will turn off automatically after approx. 5 seconds.

 **WARNING:** If the display is not on «**0.0**» ⑧, it will not be possible to figure out your weight and calculate the fat and water percentage. If the «**ERR**» ⑲ sign turns on, it means that an error has taken place during the measurement.  
Start the entire procedure again from the very beginning.

## 6. Safety, Care, Accuracy Test and Disposal

### Safety and protection

- This instrument may be used only for the purpose described in this booklet. The manufacturer cannot be held liable for damage caused by incorrect application.
- This instrument comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the «Technical Specifications» section!
- Protect it from:
  - water and moisture
  - extreme temperatures
  - impact and dropping
  - contamination and dust
  - direct sunlight
  - heat and cold
- Do not use the instrument if you think it is damaged or notice anything unusual.
- Never open the instrument.
- If the instrument is not going to be used for a prolonged period the battery should be removed.
- Read the additional safety instructions in the individual sections of this booklet.

 Ensure that children do not use the instrument unsupervised; some parts are small enough to be swallowed.

### Instrument care

Clean the instrument only with a soft, dry cloth.

### Disposal

 Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

## 7. Guarantee

This instrument is covered by a **2 year guarantee** from the date of purchase. The guarantee is valid only on presentation of the guarantee card completed by the dealer (see back) confirming date of purchase or the receipt.

- The guarantee covers only the instrument. The battery is not included.
- Opening or altering the instrument invalidates the guarantee.
- The guarantee does not cover damage caused by improper handling, a discharged battery, accidents or non-compliance with the operating instructions.

Please contact Microlife-Service (see foreword).

## 8. Technical Specifications

- Capacity: 150 kg / 330 lb / 24 st
- Resolution: 0.1 kg / 0.2 lb / 1/4 st
- 10-person memory
- One 3V lithium battery, CR 2032
- Precision and repeatability with a tolerance range within:  
+/- 1 % + 0.1 kg / 0.2 lb / 1/4 st

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